



## High-throughput DNA purification from Oragene/saliva samples<sup>†</sup> with the Magtration<sup>®</sup> 12GC

P. Lem<sup>‡</sup>, J. Chartier<sup>‡</sup>, J. Saini<sup>\*\*</sup>, and K. Obata<sup>††</sup>

<sup>‡</sup> DNA Genotek, Ottawa, Ontario, Canada

<sup>\*\*</sup> Cortex Biochem Inc., San Leandro, CA

<sup>††</sup> PSS Bio Instruments Inc., Pleasanton, CA, USA

2006-06-29

*The Magtration<sup>®</sup>12GC is an automated high-throughput DNA purification system that uses paramagnetic-particle technology. The 12GC can purify up to twelve Oragene<sup>®</sup>/saliva samples in 30 minutes, with excellent DNA yields.*

### Introduction

Large-scale population studies may involve the collection of thousands of patient samples. Manual purification of DNA from these samples can be time- and labor-intensive. The Magtration 12GC from PSS Bio Instruments is an automated DNA purification robot that uses pre-filled MagaZorb<sup>®</sup> cartridges for fast and convenient DNA purification. The purpose of this study was to determine the DNA yield of saliva samples collected with Oragene and processed with the 12GC.

### Materials and methods

#### DNA collection

Oragene self-collection kits (containing 2 mL of Oragene solution) were used to collect 2 mL of saliva from 20 donors. Prior to purification with the 12GC, the Oragene/saliva samples were incubated overnight at 50°C.

#### Automated DNA purification

DNA was extracted from 200 µL of each Oragene/saliva sample using MagaZorb DNA common kit-200 cartridges on the 12GC system. An integrated circuit


card with the DNA Common-200 protocol was supplied with the instrument. The elution volume was 200 µL.

#### DNA analysis

Purified DNA was quantified by absorbance at 260 nm. The  $A_{260/280}$  ratio was also determined.

### Results

Figure 1 shows DNA yields for the 20 Oragene/saliva samples. The median DNA yield was 3.8 µg per 200 µL of starting sample and the median  $A_{260/280}$  ratio was 1.95.



**Figure 1:** Scattergram of DNA yields from 200 µL of Oragene/saliva sample. The horizontal line represents the median yield — 3.8 µg.

<sup>†</sup> Saliva samples were collected with Oragene<sup>®</sup>-DNA or Oragene<sup>®</sup>-DISCOVER.

