

Specification: Capillary electrophoresis (CE) chip (1)

Product Code	CH_T3550_PACK	CH_T8050_PACK
Product Name	T3550	T8050
Number of chips per pack	3	2
Distance between channel and top surface	1100 μm	
Distance between channel and bottom surface	700 μm	
Channel location	No	
Total chip thickness	1800 μm	
Chip size	-	-
Channel width	50 μm	
Channel height	20 μm	
Internal Volume	-	
Number of Inlets	3	
Number of outlets	1	
Inlet/outlet hole sizes on top of the chip	1.70 mm	
Inlet/outlet holes size at channel	0.60 mm	
Optical properties	Optical clear view from all sides	
Supplied in Fluidic slide?	No	
Material chip	Borosilicate glass	



Specification: Capillary electrophoresis (CE) chip (2)

Product Code	CH_X3550_PACK	CH_X8050_PACK
Product Name	X3550	X8050
Number of chips per pack	3	2
Distance between channel and top surface	1080 μm	
Distance between channel and bottom surface	700 μm	
Channel location	No	
Total chip thickness	1800 μm	
Chip size	45 mm x 15 mm	90 mm x 15 mm
Channel width	50 μm	
Channel height	20 μm	
Internal Volume	-	
Number of Inlets	3	
Number of outlets	1	
Inlet/outlet hole sizes on top of the chip	2.50 mm	
Inlet/outlet holes size at channel	1.45 mm	
Optical properties	Optical clear view from all sides	
Supplied in Fluidic slide?	No	
Material chip	Borosilicate glass	

