

Chip Holder

Product Code	FC_PRO_CH4515	FC_PRO_CH	FC_FC4515	MR_FRAME_1507
Dimensions	128 mm x 85,4 mm x 20 mm		80 mm x 55 mm x 9.5 mm	150 mm x 75 mm x 20 mm
Material	Aluminium		Stainless Steel	
Sealing mechanism	Load 'n Seal		Slide 'n Connect	-
Maximum operation temperature	80°C		80°C (Teflon), 50°C (Fused Silica)	-
Max. operating pressure	10 bar		100 bar	10 bar, (145 psi)
Sealing material	Perlast (FFKM)		PEEK, FEP	Perlast (FFKM)
Chip thickness to fit in this holder	1.8mm excluding the black cartridge	Depending on your insert	1.8mm excluding the black cartridge	4.4mm

Product Code	746	750	883
Dimensions	128 mm x 85,4 mm x 20 mm		
Material	Aluminium		
Sealing mechanism	Load 'n Seal		
Maximum operation temperature	80°C		
Max. operating pressure	10 bar		
Sealing material	Perlast (FFKM)		
Chip thickness to fit in this holder	2.5mm after compression		

