

Specification: Cross Channel Chip

Product Code	FC_X3550CH.2_PACK	FC_X3550CH.3_PACK
Number of chips per pack	3	
Distance between channel and top surface	1100 μm	1080 μm
Distance between channel and bottom surface	680 μm	145 μm
Channel location	No	
Total chip thickness	1800 μm	1245 μm
Chip size	45 mm x 15 mm	
Channel width	50 μm	
Channel height	20 μm	
Internal Volume	0.1 μl	
Number of Inlets	3	
Number of outlets	1	
Inlet/outlet hole sizes on top of the chip	1.70 mm	
Inlet/outlet holes size at channel	0.60 mm	
Optical properties	Optical clear view from all sides	
Supplied in Fluidic slide?	Yes	
Material chip	Borosilicate glass	
Material black cartridge	Polypropylene	

