

Specification: Micromixers

Product Code	FC_SW11_PACK	FC_TD26_PACK
Number of chips per pack	2	
Distance between channel and top surface	950 μm	
Distance between channel and bottom surface	550 μm	
Channel location	No	
Total chip thickness	1975 μm	1800 μm
Chip size	45 mm x 15 mm	
Channel width	200 μm	
Channel height	150 μm	
Internal Volume	1 μl	2 μl
Number of Inlets	2	
Number of outlets	1	
Inlet/outlet hole sizes on top of the chip	1.70 mm	
Inlet/outlet holes size at channel	0.60 mm	
Optical properties	Clear sight from the top of chip only	
Supplied in Fluidic slide?	Yes	
Material chip	Borosilicate glass	
Material black cartridge	Polypropylene	

